

400 MHz Photoreceiver with Si-PIN Photodiode



The picture shows model HCA-S-400M-SI-FS. The photoreceiver will be delivered without post holder and post.


Features

- Si-PIN photodiode, 0.8 mm active diameter
- Bandwidth DC 400 MHz
- Amplifier transimpedance gain 5.0 × 10³ V/A
- Max. conversion gain 2.7 x 10³ V/W @ 800 nm
- Spectral range 320 1000 nm
- Free-space input 1.035"-40 threaded, alternatively 25 mm diameter unthreaded
- Easily convertible to fiber optic input (FC and FSMA) with optionally available screw-on adapters
- Fiber optic input also available as permanently mounted FC- or FSMA-input
- UNC 8-32 and M4 tapped holes for mounting on standard posts with metric and imperial thread

Applications

- Spectroscopy
- Fast pulse and transient measurements
- Optical triggering
- Optical front-end for oscilloscopes, A/D converters and HF lock-in amplifiers

Block Diagram



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O

HCA-S-400M-SI_R5/TH,JMa/04MAY2022 Page 1 of 11

Available Versions

HCA-S-400M-SI-FST



Picture shows 1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm)

1.035"-40 threaded flange for free space applications. Compatible with many optical standard accessories and for use with various types of fiber connector adapters.

Optionally available:

Fiber adapters PRA-FC, PRA-FCA and PRA-FSMA. The coupling efficiency will depend on fiber type. With the relative large 0.8 mm dia. photodiode installed in the HCA-S-400M-SI input coupling is not critical. However, standard SM 9/125 fibers (PC or APC) with low numerical aperture (NA) are recommended for ensuring near 100% coupling efficiency.

HCA-S-400M-SI-FS



Picture shows unthreaded flange with 25 mm diameter

25 mm dia. unthreaded flange for free space applications. Compatible with many optical standard accessories.

HCA-S-400M-SI-FC



Fix/permanent FC fiber connector for high coupling efficiency and excellent conversion gain accuracy.

HCA-S-400M-SI-SMA



Fix/permanent FSMA fiber connector for high coupling efficiency and excellent conversion gain accuracy.

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E TO

400 MHz Photoreceiver with Si-PIN Photodiode

Related Models HCA-S-400M-IN-FST InGaAs-PIN, \varnothing 0.3 mm, 900 - 1700 nm free space input, 1.035"-40 threaded flange

HCA-S-400M-IN-FS InGaAs-PIN, Ø 0.3 mm, 900 - 1700 nm free space input, 25 mm dia. unthreaded flange

HCA-S-400M-IN-FC InGaAs-PIN, integrated ball lens, 900 - 1700 nm

FC fiber connector (fix/permanent)

Available Accessories PRA-FC

PRA-FCA PRA-FSMA







Fiber-adapter with external 1.035"-40 thread (suitable for FST models only).

PRA-PAP



Alternative mounting option: Post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, PWPR, HCA-S and LCA-S.

PS-15-25-L



Power supply Input: 100 – 240 VAC Output: ±15 VDC

Specifications Test conditions

 $V_S=\pm 15$ V, $T_A=25$ °C, output load impedance 50 Ω , warm-up 20 minutes (min. 10 minutes recommended)

Gain Transimpedance gain

gain 5.0 \times 10³ V/A (@ output load 50 Ω) \pm 1 % (electrical)

Gain accuracy Conversion gain

 2.7×10^3 V/W typ. (@ 800 nm, output load 50 Ω)

Frequency Response

Lower cut-off frequency DC

Upper cut-off frequency (–3 dB) 400 MHz (±10 %)

Gain flatness

±1 dB

Time Response

Rise/fall time (10 % – 90 %)

1.0 ns

Input

Noise equivalent power (NEP) Optical saturation power Input offset compensation range 40 pW/ $\sqrt{\text{Hz}}$ (@ 800 nm, 100 MHz) 400 μ W (for linear amplification, @ 800 nm) \pm 200 μ A, adjustable by offset potentiometer

Detector

 $\begin{array}{lll} \mbox{Detector} & \mbox{Si-PIN photodiode} \\ \mbox{Active area} & \mbox{\varnothing 0.8 mm} \\ \mbox{Spectral range} & \mbox{320} - 1000 \mbox{ nm} \\ \end{array}$

Max. sensitivity

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

0.55 A/W typ. (@ 800 nm)

F E T O

Datasheet

HCA-S-400M-SI

400 MHz Photoreceiver with Si-PIN Photodiode

Specifications (continued)			
Output	Output voltage range Max. output voltage range Output impedance Output noise	± 1.0 V (@ 50 Ω output for linear operation and ± 1.5 V (@ 50 Ω load) 50 Ω (terminate with 50 3 mV _{RMS} (20 mV _{PP}) typ. (detector, measurement	low harmonic distortion Ω load) Ω load, no signal on
Input Flange	Material	1.4305 stainless steel, nickel-plated (FST flange) AlMg4.5Mn, nickel-plated (FS flange)	
Coupler Ring (FST version only)	Material	1.4305 stainless steel, (glass bead blasted
Power Supply	Supply voltage Supply current	± 15 V (± 14.5 V ± 16.5 V) ± 55 mA (depends on operating conditions, recommended power supply capability min. ± 150 mA)	
Case	Weight Material	209 g (0.46 lbs) HCA-S-400M-SI-FST incl. coupler ring 196 g (0.43 lbs) HCA-S-400M-SI-FS 188 g (0.41 lbs) HCA-S-400M-SI-FC 200 g (0.44 lbs) HCA-S-400M-SI-SMA AlMg4.5Mn, nickel-plated	
Temperature Range	Storage temperature Operating temperature	-30 °C +85 °C 0 °C +60 °C	
Absolute Maximum Ratings	Optical input power (CW) Power supply voltage	20 mW ±20 V	
Connectors	Input	HCA-S-400M-SI-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories
		HCA-S-400M-SI-FS	25 mm dia. unthreaded flange for free space applications
		HCA-S-400M-SI-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible)
		HCA-S-400M-SI-SMA	FSMA fiber optic connector (fix/permanent)
	Output	BNC jack (female)	
	Power supply	LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)	
		PIN 2 O O O	PIN 1 +Vs Pin 1: +15 V Pin 2: -15 V PIN 3 GND Pin 3: GND
Scope of Delivery	HCA-S-400M-SI, internally threaded coupler ring (FST version only), LEMO® 3-pin connector, datasheet, transport package		


SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E T O

Datasheet HCA-S-400M-SI **400 MHz Photoreceiver** with Si-PIN Photodiode Ordering Information HCA-S-400M-SI-FST 1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories. HCA-S-400M-SI-FS 25 mm dia. unthreaded flange for free space applications. HCA-S-400M-SI-FC FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible). HCA-S-400M-SI-SMA FSMA fiber optic connector (fix/permanent). Spectral Responsivity 0.6 0.5 Photo Sensitivity (A/W) 0.4 0.3 0.2 0.1 0 600 700 800 200 300 400 500 900 1100 1000 Wavelength (nm) DB-Sens-HCA-S-400M-SI_R01 Typical Performance Frequency response Characteristics RBW 3 MHz Att 0 dB * VBW 1 MHz -3.36 dB M1[1] Ref -20.0 dBm SWT 2.5ms 408.440000000 MHz 10 dB 5 dB 0 dB -5 dB--10 dB -15 dB -20 dB -25 dB--30 dB Start 20.0 MHz Stop 800.0 MHz PD-HCA-S-400M-Si-bw_R01 SOPHISTICATED TOOLS FOR SIGNAL RECOVERY 0

HCA-S-400M-SI_R5/TH,JMa/04MAY2022 Page 5 of 11

Typical Performance Characteristics (continued)




Note: spectral noise data is measured at the amplifier output with no signal on the photodiode. To determine the spectral input noise divide the measured output noise by the amplifier conversion gain.

Conversion gain (V/W) = amplifier gain (V/A) \times photo sensitivity (A/W).

Marker	frequency	output noise	resulting input noise (NEP)
1	100 MHz	109 nV/√Hz	40 pW/√Hz (@ 800 nm)

Pulse response to square wave input signal (with 16 times averaging)




PD-HCA-S-400M-Si pulse-2ns_R01

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E TO


400 MHz Photoreceiver with Si-PIN Photodiode

Typical Performance Characteristics (continued) Large signal response output signal for 100 MHz, 370 μ W modulated optical input signal (with 4 times averaging)



PD-HCA-S-400M-Si-large-sinus_R01

Small signal response output signal for $3.7~\mu W$ modulated optical input signal, 1 MHz square wave, without averaging



PD-HCA-S-400M-Si noise-square_R01

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY


F E M T O

Datasheet

HCA-S-400M-SI


400 MHz Photoreceiver with Si-PIN Photodiode

Typical Performance Characteristics (continued) Small signal response output signal for 3.7 μ W modulated optical input signal, 1 MHz square wave, with 64 times averaging



Dimensions

HCA-S-400M-SI-FST (1.035"-40 threaded free space input)





SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O

Dimensions (continued)

HCA-S-400M-SI-FS (25 mm dia. unthreaded free space input)





DZ-HCA-S-X00-SI_FS_R1


all dimensions in mm unless otherwise noted

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E MT O

Dimensions (continued)

HCA-S-400M-SI-FC (FC fiber optic connector)





DZ-HCA-S-XX-XX_FC_R1

all dimensions in mm unless otherwise noted

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O







400 MHz Photoreceiver with Si-PIN Photodiode

Dimensions (continued)

HCA-S-400M-SI-SMA (FSMA fiber optic connector)





DZ-HCA-S-XX-XX_SMA_R1

all dimensions in mm unless otherwise noted

FEMTO Messtechnik GmbH Klosterstr. 64 10179 Berlin · Germany Phone: +49 30 280 4711-0 Fax: +49 30 280 4711-11 Email: info@femto.de www.femto.de Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O

HCA-S-400M-SI_R5/TH,JMa/04MAY2022 Page 11 of 11